SORBACT® COMPRESS

EN Instructions for Use

SV Bruksanvisning

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DA Brugsanvisning

NO Bruksanvisning

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ET Kasutusjuhised

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AU/NZ Bayport Brands Pty Ltd, 21 Bayport Court, Mornington. Vic 3931, Australia, Ph: +61 59 25 90 00, sales@bayportbrands.com.au, www.bayportbrands.com.au

TR NIBYKAN Golden Pharma Mannheim GmbH, Augustaanlage 32, 68165, Mannheim, Germany, Tel. +49 621 43 17 93 39

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17 cm x 28 cm
7 in x 11 in

11 cm x 15 cm
4 in x 6 in

Read the instructions for use before starting treatment with Sorbact® Compress.

Device Description

Sorbact® Compress is a non-absorbent, bacteria and fungi binding wound contact layer, based on Sorbact® technology, that allows passage of wound exudate into a secondary dressing. Sorbact® Compress prevents and treats wound infections and facilitates the wound healing process.

Intended Purpose

Sorbact® Compress is intended for use in management of clean, contaminated, colonized or infected exuding wounds, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic foot ulcers and leg ulcers. Sorbact® Compress can be used on both superficial and deep wounds.

Intended User and Use Environment

Sorbact® Compress is intended to be used on children and adults. The dressing can be used by healthcare professionals, or by lay persons under supervision of a healthcare professional. The dressing is intended to be used in healthcare facilities and home environment.

Mode of Action

Sorbact® Compress binds common wound microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus* species, *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. There is no release of antimicrobial agents to the wound.

Clinical Benefit

Use of Sorbact® Compress can reduce bioburden and prevent infection. The dressing helps reducing odor and pain, and can improve wound healing and decrease wound size.

Contraindications

There are no known contraindications to the use of Sorbact® Compress. Warnings and precautions should be noted.

Warnings and Precautions

Do not use on patients with a known hypersensitivity to the dressing.

Do not use in combination with fatty products, such as ointments, creams and solutions as they may decrease the binding of microorganisms.

Do not reuse. Sorbact® Compress is for single patient and single use only. Reuse may lead to cross contamination and increased risk of infection.

Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed.

Do not wrap around extremities, as circulation may be compromised.

Do not apply Sorbact® Compress in direct contact with bone and tendons, exposed blood vessels, organs or nerves.

Do not re-sterilize.

If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

Instructions for Use

1. Prepare the wound and surrounding skin according to local clinical practice.
2. Select an appropriate dressing size for the wound. The dressing can be unfolded into suitable size and may overlap the wound margins if needed.
3. Remove the dressing from the peel pouch using an aseptic technique.
4. If the dressing is cut, use an aseptic technique. Do not tear the dressing. Discard any unused dressing.
5. Apply the dressing. Ensure that the dressing comes into direct contact with the complete wound surface to allow microorganisms to bind to the dressing. In deep wounds, unfold and fluff up one or more dressings and fill the wound. Avoid dense packing. The dressing can also be used as a liner.
6. Apply a secondary dressing appropriate for the exudate level and fixate.
7. The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.
8. Should the dressing adhere to the wound, moisten the dressing to assist removal and to avoid disruption of the healing wound.

Sorbact® Compress is MRI safe.

Sorbact® Compress can be used in conjunction with compression therapy.

Sorbact® Compress can be used in conjunction with NPWT. Follow the Instructions for Use of the NPWT system.

Storage and Disposal

Sorbact® Compress shall be stored dry and kept away from sunlight.

Disposal should be made according to local environmental procedures.

Notice to User

Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB, (complaints@abigo.com) and the competent authority of your state.

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