PRODUCT SPECIFICATION

REF 51011005 BURN GEL 100ML

1 GENERAL PRODUCT DESCRIPTION

Description

Cederroth Burn Gel 100ml is a water-based gel for minor burns, scalds, sunburn and radiation damage. It acts as a temporary protective barrier that cools and soothes the skin and thereby provides pain relief.

The gel is simple to apply (gel in an aerosol dispenser) and it can be used multiple times, as it is protected against air until last drop.

Instructions for use: Apply the gel generously on the affected area as soon as possible after the injury. Can also be used on children. Easy to wash off.

Precautions: Do not use on deeper second degree burns or other types of burns such as chemical burns, phototoxic reactions, skin induced drug reactions or skin damages induced by electrical chock.

Contains preservatives. CE labelled medical device, class IIa.

Storage information

Dry storage in room temperature. Do not expose to temperatures exceeding 50°C/122°F. Also refer to warnings in labelling.

2 DETAILED PRODUCT INFORMATION

Content per consumer unit

<u>Product/Material</u> Burn Gel	!	<u>Quantity</u> 100 ml	<u>Dimensions (d x h mm)</u> Can incl. lid: 45 x 154
Material		100 111	
	Purified water		Sodium Propyl parabon
Liquid:	Pullieu water		Sodium Propyl paraben
	Carbomer 980NF		Aloe Vera (200:1 dry powder)
	Sodium Methyl paraben		NaOH (20% w/w in water)
Propellant:	Nitrogen		
Bag-on-valve:	Aluminum bag laminated wit	th PP-OPA-ALU-PET	
Actuator:	PP (Polypropylene)		
Cap:	PP (Polypropylene)		
Can:	Aluminum		
Multi-label:	Polypropylene, white gloss fi	lm.	

3 LOGISTIC INFORMATION

Shelf life:

	CU	SKU	Pallet	
GTIN:	7310610015275	7310800015436	7312020016040	
Content:	1 aerosol	8 CU	48 SKU/layer x 6 layers	
Size (w x l x h mm):	45 x 45 x 154	191 x 97 x 171	800 x 1200 x 1176	
Gross weight:	135 g	1145 g	355 kg	
Packing material:	-	Corrugated board 65 g	-	
Transport information:		UN1950, non-flammable a Limited quantity LQ	UN1950, non-flammable aerosols, class 2.2, Limited guantity LQ	



Version: 01 Issued: 2019-01-11





Orkla Care AB, Svetsarvägen 15, Box 1336, SE 171 26 Solna, Sweden, Phone +46 (0)10 142 64 00 www.cederroth.com firstaid@cederroth.com

3 years from production date